




# Migration Patterns of DNA Ladder Markers in Different % of Agarose Gels



1. PCR Marker 0.1-1Kbp (EBM-1005)
2. PCR Marker 0.1-3Kbp (EBM-1006)
3. PCR Marker plus 50bp (EBM-1007)

0.7%, 1%, 1.5%  
 0.5x TBE Gel  
 100V constant  
 EtBr staining



1. 1Kbp DNA Mass Marker (EBM-1009)
2. 100bp DNA Mass Ladder Marker (EBM-1010)
3. 500bp Step Ladder Marker (EBM-1101)
4. 200bp Step Ladder Marker (EBM-1102)
5. 100bp Step Ladder Marker (EBM-1103)
6. 50bp Step Ladder Marker (EBM-1104)

1.5%, 0.5x TBE Gel  
 100V constant  
 EtBr staining